

## Method: Face, Ear and Nec

## EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	1Hz	10	60	1.2
	Brown	1Hz	12	70	1.2
	Blonde	1Hz	14	90	1.2
3&4	Dark Brown/Black	1Hz	8	60	1.2
	Brown	1Hz	10	80	1.2
5&6	Dark Brown/Black	1Hz	6	100	1.2
	Brown	1Hz	6	100	1.2

## Method: Abdomen, Chest

## EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	4Hz	9	60	10
	Brown	4Hz	10	70	10
	Blonde	4Hz	11	90	10
3&4	Dark Brown/Black	4Hz	9	60	10
	Brown	4Hz	10	80	10
5&6	Dark Brown/Black	4Hz	4	100	10
	Brown	4Hz	6	100	10

## Method: Underarm

## EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	4Hz	9	60	3.2
	Brown	4Hz	10	70	3.2
	Blonde	4Hz	11	90	3.2
3&4	Dark Brown/Black	4Hz	9	60	3.2
	Brown	4Hz	10	80	3.2
5&6	Dark Brown/Black	4Hz	4	100	3.2
	Brown	4Hz	6	100	3.2

## Method: Arm and Hands

## EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	4Hz	9	60	9.6
	Brown	4Hz	10	70	9.6
	Blonde	4Hz	11	90	9.6
3&4	Dark Brown/Black	4Hz	9	60	9.6
	Brown	4Hz	10	80	9.6
5&6	Dark Brown/Black	4Hz	4	100	9.6
	Brown	4Hz	6	100	9.6

Method: Inside crease at top of thigh

EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	1Hz	9	60	3.2
	Brown	1Hz	11	70	3.2
	Blonde	1Hz	13	80	3.2
3&4	Dark Brown/Black	1Hz	9	60	3.2
	Brown	1Hz	11	80	3.2
5&6	Dark Brown/Black	1Hz	4	100	3.2
	Brown	1Hz	6	100	3.2

Method: Brazilian / Hollywood

EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	1Hz	9	60	3.2
	Brown	1Hz	11	70	3.2
	Blonde	1Hz	13	80	3.2
3&4	Dark Brown/Black	1Hz	9	60	3.2
	Brown	1Hz	11	80	3.2
5&6	Dark Brown/Black	1Hz	4	100	3.2
	Brown	1Hz	6	100	3.2

Method: Leg, Knee and Foot

EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	4Hz	9	60	7.5
	Brown	4Hz	10	70	7.5
	Blonde	4Hz	11	90	7.5
3&4	Dark Brown/Black	4Hz	9	60	7.5
	Brown	4Hz	10	80	7.5
5&6	Dark Brown/Black	4Hz	4	100	7.5
	Brown	4Hz	6	100	7.5

Method: Back and Nape

EVO VELOCITY TRIO

Skin Type	Natural Hair colour	Frequency	Energy Setting	Pulse Width	Total KJ/Area
1&2	Dark Brown/Black	4Hz	9	60	10
	Brown	4Hz	10	70	10
	Blonde	4Hz	11	90	10
3&4	Dark Brown/Black	4Hz	9	60	10
	Brown	4Hz	10	80	10
5&6	Dark Brown/Black	4Hz	4	100	10
	Brown	4Hz	6	100	10